


COLOR TV CAMERA FOR MICROSCOPES

CK3900/3100

INSTRUCTION MANUAL



MEIJI TECHNO CO., LTD.

JAPAN



DECLARATION OF CONFORMITY

AS DEFINED BY THE COUNCIL DIRECTIVE

89/336/EEC
EMC (ELECTROMAGNETIC COMPATIBILITY)

WE HEREWITH DECLARE THAT THIS PRODUCT
COMPLIES WITH THE FOLLOWING PROVISIONS APPLYING TO IT.

EN-50081-1

EN-50082-1

ENGLISH VERSION



The lightning flash with arrowhead symbol, within an equilateral triangle, is intended to alert the user to the presence of uninsulated "dangerous voltage" within the product's enclosure that may be of sufficient magnitude to constitute a risk of electric shock to persons.



The exclamation point within an equilateral triangle is intended to alert the user to the presence of important operating and maintenance (servicing) instructions in the literature accompanying the appliance.



The slash within a circle is intended to alert the user to the presence of prohibition of any kind of operation, maintenance and storage.



The abstraction within a circle is intended to alert the user to presence of prohibition to disassemble the product.



The abstraction within a circle is intended to alert the user to turn the power off on the product and to take out the plug.



The abstraction within a circle intended alert the user to presence of prohibition to expose the product to rain, moisture or any kind of moist surroundings.

Precautions



Do not attempt to disassemble this camera.



To prevent electric shock, do not remove cover. There are no user-serviceable parts inside. Refer servicing to qualified service personnel.



Do not expose this camera to rain or moisture.



Do not face this camera towards the sun, extreme bright light or light reflecting objects.

Even when this camera is not in use, put the supplied lens cap on the lens mount.



Handle this camera with the maximum care.



Operate this camera only from the type of power source indicated on the camera.



Power-off the camera during any modification such as changes of jumper-line and jumper-register.

USER'S RECORD

The production serial number are shown on the bottom of camera

Model Name :

Serial No. :

Table of Contents

CE declaration.....	2
Caution.....	3
Precaution.....	4
1. General.....	6
2. Main features.....	6
3. Standard composition.....	7
4. Pin assignment.....	7
5. Functions and operations.....	8
5.1 SW1 switch on PK8167 board.....	8
5.2 On-screen-menu.....	9
5.3 RS 232C functio ns ns.....	13
6. External appearance.....	14
7. Specifications.....	15
Typical CCD characteristic.....	16

1. General

The CV-S3200/3300 is a new high-resolution color camera series, available in PAL or NTSC format. The camera series was developed with JAI's new DSP processor, which has remarkable high signal to noise ratio as well as a high bit flow. The cameras also have built-in RAM, which together with the DSP enables the camera to produce field integration with 100% horizontal resolution.

The CV-S3200/3300 series features EXview HAD™ sensors (CV-S3200:1/2" and CV-S3300:1/3") which are extra light sensitive and have enhanced anti-smear and blooming characteristics.

The cameras run on 12V DC and can be gen-locked or use VBS for external synchronization. There are connections for video control or direct drive lens, and VBS and Y/C video outputs. Both cameras have CS/C-mount with a variable mounting ring for back focus. The CV-S3300 has screw terminals for power, gen-lock and line-lock. The CV-S3200 comes with D-sub connector.

The camera set-up can be performed via the JAI software platform and an RS232C port, or from a new on-screen menu that can be controlled by the three buttons mounted on the rear of the camera. Both the software and the on-screen-menu can control the following setting : shutter, integration, gamma, AGC, H/V enhancer, iris-level for direct drive, positive/negative picture, back-light compensation, and AWB fixed or manual color temperature.

Applications include surveillance, security, medical and a large variety of other industrial applications.

2. Main Features

- 1/2" Exview HAD™ IT color CCD sensor for CV-S3200, 1/3" Exview HAD™ IT color CCD sensor for CV-S3300
- NTSC : 768 (h) x 494 (v) , PAL : 752 (h) x 582 (v)
- 450 TV lines horizontal resolution
- Composite VBS and Y/C output
- Excellent color reproduction
- Frame integration mode for higher sensitivity
- S/N ratio > 50 dB
- White balance automatic or manual
- Auto iris and CCD iris for light regulation
- Manual shutter up to 1/10,000 sec.
- Scanning reversal function (mirror-effect mode)
- CS-lens mount with focus adjustment
- CV-S3200 employs 9 pin Dsub connector at camera rear panel for single cable solution.
- CV-S3300 provides 7 pin terminal at camera rear panel for easy cable installation.


3. Standard compositions

- | | | |
|----|-----------------------------------|-----|
| 1) | Camera main body | x 1 |
| 2) | C-mount adaptor (MP-55) | x 1 |
| 3) | Auto iris lens pulg (E4-191)-100) | x 1 |
| 4) | Y/C out plug (TCP8040-01-520) | x 1 |
| 5) | Operation manual | x 1 |

4. Pin assignment


4-1. DC IN/SYNC/RS-232C connector

a) 9 pin D-sub connector for CV-S3200




- 1 : Ground
- 2 : +12V DC input
- 3 : Ground
- 4 : VBS output
- 5 : Sync. input
- 6 : Y output
- 7 : C output
- 8 : TXD output
- 9 : RXD input

b) 7 pin terminal for CV-S3300




- 1 : RXD input
- 2 : TXD output
- 3 : Ground
- 4 : Sync. input
- 5 : Ground
- 6 : +12V DC input
- 7 : Ground

4-2. S-VIDEO connector



- 1 : Ground
- 2 : Ground
- 3 : Y output
- 4 : C output

4-3. LENS connector




- 1 : +9V / Damp-
- 2 : NC / Damp+
- 3 : Iris / Drive+
- 4 : Ground / Drive-

5. Functions and operations

5-1. SW1 switch on PK8167 board

This switch is to select the type of auto iris. It is located on the CCD board (PK8167), as below.

- Up : Direct drive mode (Passive)
- Down : Video controlled (Active)



Note : Factory pre-set is Video controlled.

5-2. On-screen-menu

On-screen-menu can be controled by the "UP", "DOWN" and "FUNC" switches on the rear panel.

5-2-1. Functions

The following functions are adjustable at on-screen-menu.

- White balance
- Shutter
- Gain
- Back light compensation
- Lens iris
- Gamma
- Scanning reversal
- Sensitive

a) White balance

This function is to adjust the color temperature to make the color of the image correct.

3200 K : For incandescent lamp or halogen lamp.

4600 K : For white fluorescent lamp.

5600 K : For daylight.

Auto : The TTL system is employed for auto white balance, and it detects the highest video level in objects to recognize a white color.
It is necessary to place the camera in a position where the camera does not pick up spot lights and fluorescentlights.

Manual : To adjust the color temperature between 2700 K to 6700 K.

Note : Factory pre-set is auto white balance.

b) Shutter

This function is used for "freezing" objects and for light regulations. Nine exposure time settings are available.

OFF : 1/60 sec. (NTSC), 1/50 sec. (PAL)

Selectable shutter : 1/125, 1/250, 1/500, 1/1000, 1/2000, 1/4000, 1/10,000 sec.

Flicker-less : 1/100 sec. (NTSC), 1/120 sec. (PAL)

CCD iris : Shutter speed changes automatically.
(1/60 sec. to 1/10000 sec. (NTSC), 1/50 sec. to 1/10,000 sec. (PAL))

Note : Factory pre-set is OFF.

CAUTION

1. Under certain unique light conditions, the following may appear.
 - Strong smear and/or blooming on highlighted objects.
 - Noticeable flicker in the picture.
2. CCD iris
Continuous variable electronic shutter controls the exposure time of the CCD sensor, depending on the incoming light level. When this mode is selected, auto or manual iris lens can be used instead of an auto iris lens.
3. Shutter speed and sensitive mode
In sensitive mode the CCD iris does not function.

c) AGC (Auto gain control)

OFF : Gain is fixed at minimum 0 dB

+12 dB : AGC ON (0 to +12 dB)

+24 dB : AGC ON (0 to +24 dB)

Manual : Manual gain

Note : Factory pre-set is AGC OFF.

CAUTION

1. Auto iris lens can be used only when the CCD iris is OFF.

d) BLC (Back Light Compensation)

BLC circuite **will** compensate for extreme highlights in the field of view.

There are **three** pre-defined patterns that can be selected. These are shown below. The higher the value, the **more** backlight compensation.

1	1	
4	9	4
1	1	

0	0	
0	15	0
0	0	

0	0	
0	10	0
4	4	

Note : Factory pre-set is OFF.

CAUTION

- The BLC **works** under the follwing conditions:

BLC	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
AGC	ON	ON	OFF	OFF
CCD iris	ON	OFF	ON	OFF

: Effective
 : Not-effective

- It is **recommended** to switch "BLC OFF" when using the camera with pan/tilt head or if the scene illumination changes **rapidly**.
- With "BLC **ON**" the regulation speed of the iris is decreased.

e) Iris level

This is to set **the** video level, when the camera is set at CCD iris, AGC ON and auto iris.

f) Gamma

This function is to set the gamma correction.

- 0.45 : Correction
- 0.6 : Correction
- 1.0 : Linear

Note : Factory pre-set is 0.45.

g) Scanning reversal (mirror-effective mode)

This function can change the horizontal scanning between normal and reverse.

- OFF : Normal scanning
- ON : Reversal scanning

Note : Factory pre-set is set at OFF.

h) Frame integration


The sensitivity can be set at these pre-set positions.

- OFF : Standard sensitivity (Exposure time = 1 field)
- Sensitive 1 : 2 times higher than standard sensitivity (Expoure time = 2 fields)
- Sensitive 2 : 4 times higher than standard sensitivity (Exposure time = 4 fields)


Note : Factory pre-set is OFF.

5-2-2. Operations

Screen display (example)



Function buttons for on-screen-menu



Rear panel of CV-S3300 (example)

Note : This scroll bar will be displayed when the following functions are selected.

White balance = manual, AGC = manual, Iris = manual

- (a) To start adjustment, push the "FUNC" button. The menu will appear on the monitor .
(See page 12)
- (b) Use "UP" or "DOWN" buttons to choose functions. To select function, press the "FUNC" button.
- (c) After selecting the function, press "UP" or "DOWN" buttons to choose mode. To select mode, press the "FUNC" button.
- (d) Go to the next function if needed, and repeat (b) to (c).


5-3. RS 232C functions

- a) White balance : To adjust the color temperature
- b) Electronic shutter : To adjust the shutter speed
- c) Auto gain control : OFF, +12dB, +24dB, manual selectable
- d) Back light compensation : To create your own BLC pattern.
- e) Iris level : To adjust the video level
- f) Gamma : 0.45, 0.60, 1.0 selectable
- g) Scanning reversal : To select the scanning direction
- h) Frame integration : To select the sensitivity


6. External appearance

Unit : mm (inches)

e.g. CV-S3200



e.g. CV-S3300



7. Specification

Model name	CV-S3200	CV-S3300
TV standard	NTSC / PAL	NTSC / PAL
Scanning system	NTSC : 525 lines, 30 frames/sec. PAL : 625 lines, 25 frames/sec.	
CCD sensor	Color 1/2" IT CCD sensor	Color 1/3" IT CCD sensor
Effective pixels	NTSC : 768 (h) x 494 (v) PAL : 752 (h) x 582 (v)	
Cell size	NTSC : 8.4 (h) x 9.8 (v) μm PAL : 8.6 (h) x 8.3 (v) μm	NTSC : 6.35 (h) x 7.4 (v) μm PAL : 6.3 (h) x 6.25 (v) μm
Resolution (horizontal)	450 TV line	
Sensitivity	0.4 lux on sensor (AGC ON, 50% video)	
S/N ratio	50 dB (AGC OFF, Gamma=1.0)	
Video output	VBS 1.0 Vpp at 75 Ohm Y/C output 1.0 Vpp at 75 Ohm	
Gain	Off, +12 dB, +24 dB, manual	
White balance	3200K, 4600K, 5600K, auto, manual	
Gamma	1.0 - 0.6 - 0.45	
Synchronization	Internal / external sync. or VBS 1.0 Vpp	
Shutter	OFF (1/60 at NTSC, 1/50 at PAL), 1/125, 1/250, 1/500, 1/1000, 1/2000, 1/4000, 1/10,000, CCD iris, Flicker-less (NTSC : 1/100, PAL : 1/120)	
Frame integration	Off, sensitive 1, sensitive 2	
Back light compensation	Off, pattern 1, pattern 2, pattern 3	
Scanning reversal	Off, On	
Serial interface	RS 232C	
Rs 232C functions	Electronic shutter, CCD iris, Flicker-less, White balance, Back light compensation, Video reversal, Gamma, Gain, Chroma	
Operating environment	Temperature : -5°C to +45°C Humidity : 20% to 80% non-condensing	
Power	+12 VDC \pm 10%, 3.0W	
Lens mount	CS-mount	
Dimensions	45 x 55 x 110.2 (HxWxD) mm	
Weight	350 g approx.	

Note : Above specifications are subject to change without notice.