



GEM STEREO MICROSCOPES

For the professional Jeweler/Gemologist

MEIJI TECHNO'S ADVANCED GEM STEREO MICROSCOPE SYSTEMS
FEATURE INCLINABLE, SWIVEL BASE STANDS, WITH BUILT IN
TRANSMITTED LIGHT DARKFIELD ILLUMINATION.

TRANSMITTED ILLUMINATION IS PROVIDED BY A 30 WATT HIGH INTENSITY
HALOGEN LAMP WITH A VARIABLE INTENSITY CONTOROL. INCIDENT
ILLUMINATION IS PROVIDED BY A NEW BRIGHTER
FLUORESCENT ILLUMINATION SYSTEM.



Model: GEMZ-5SVH

GEM BINOCULAR ZOOM STEREO SYSTEM

Composed of:

- EMZ-5** Binocular zoom stereo microscope body with 0.7X to 4.5X zoom range, magnification from 7X to 90X(Using SWF 10X and SWF20X eyepieces)
- MA502** Super widefield 10X eyepiece, FN 23mm(paired)
- MA504** Super widefield 20X eyepiece, FN 11.5mm(paired)
- GEMZ/SVH** Professional swivel GEM stand, inclinable and rotatable for greater ease of use, with built-in darkfield transmitted light illumination system, 30W halogen lamp with variable intensity control.
- MA557** Iris diaphragm with adapter.
- MA309** Fluorescent box illuminator with articulating arm and 115V or 220/240V transformer
- MA555/05** Gem clamp

Model: GEMT-2SVH

GEM BINOCULAR TURRET STEREO SYSTEM

Composed of:

- EMT-2SVH** Binocular stereo body with 1X and 3X objectives in a turret magnification changer. Magnification range: 10X, 20X, 30X and 60X.(Using 10X and 20X eyepieces)
- MA502** Super widefield 10X eyepiece, FN 23mm(paired)
- MA504** Super widefield 20X eyepiece, FN 11.5mm(paired)
- GEMT/SVH** Professional swivel GEM stand, inclinable and rotatable for greater ease of use, with built-in darkfield transmitted light illumination system, 30W halogen lamp with variable intensity control.
- MA557** Iris diaphragm with adapter.
- MA309** Fluorescent box illuminator with articulating arm and 115V or 220/240V transformer
- MA555/05** Gem clamp



EYEPIECES AND ACCESSORIES FOR GM SERIES


- MA502** Super widefield eyepiece 10X (paired), FN:23
- MA520** Super widefield eyepiece 12.5X (paired), FN:20
- MA503** Super widefield eyepiece 15X (paired), FN:15.4
- MA535** Super widefield eyepiece 15X high-eyepoint (paired), FN:16
- MA504** Super widefield eyepiece 20X (paired), FN:11.5
- MA521** Super widefield eyepiece 30X (paired), FN:7.7
- MA519** Super widefield focusing eyepiece 10X (ea)
- MA600** Eyeshields for SWF eyepieces (paired)

SPARE BULBS

- MA326** Halogen lamp, 6V 30W
- MA309/05** Replacement bulb for MA309/100 and MA309/200 fluorescent illuminators, 4W,
(Note: MA309/100 and MA309/200 require two MA309/05 bulbs)

REPLACEMENT COMPONENTS

- MA309/100** Fluorescent box illuminator with articulating arm for use with ABH/1 and ABH/1,115V
- MA309/200** Fluorescent box illuminator with articulating arm for use with ABH/1 and ABH/1,220/240V
- MA309/15** Replacement transformer (without external case) for MA309/100, 115V
- MA309/20** Replacement transformer (without external case) for MA309/200, 220/240V
- MA555/05** Stone holder
Spring-loaded, with dual axis rotation. For stones and mounted goods.
- MA555/10** Inclusion pointer
Moves across field of view to pinpoint gem features.



DIAMOND LOUPES

Highest quality loupe with Hasting triplet lenses fully corrected for spherical and chromatic aberration. Meets regulations for diamond grading. Supplied in an attractive leather case.

- MG600/10** Diamond loupe, 10X triplet 21.0mm diameter, agg-shaped, black anodized
- MG610/10** Diamond loupe, 10X triplet 18.0mm diameter, agg-shaped, black anodized



MG600/10

MG610/10



Agent :



MEIJI TECHNO CO.,LTD.

322-1,Chikumazawa Miyoshi-machi, Iruma-gun
Saitama 354-0043, Japan
E-mail : meiji@meijitechno.co.jp
Phone : (0)49-259-0111 Fax : (0)49-259-0113



MEIJI TECHNO AMERICA

5895 Rue Ferrari, San Jose, CA 95138, USA
Tel : 408-226-3454 Fax : 408-226-0900
Toll Free : 1-800-832-0060
E-mail : info@meijitechno.com
<http://www.meijitechno.com>