






S-4000 Series UNIVERSAL STAND

S-4000 Series Universal Stand is a sturdy stand for use in industrial work station where a wide horizontal and vertical viewing latitude is required. Due to its solid stability, no counter-balance is required.

The horizontal pole is rotatable 360° around the vertical pole and can be movable up and down along the vertical pole and fixable at any desired position by a clamp screw on the slider block. The horizontal pole is slotted to prevent accidental rotation of the Zoom Stereo Body.



S-4100 S-4200

S-4000 Series UNIVERSAL STAND

MODELS & SPECIFICATIONS

MODELS S-4100: UNIVERSAL STAND WITH POST FOR F HOLDER

BASE: 265mm(W) x 265mm(D) x 25.4mm(H)
(10.43" x 10.43" x 1.0")

VERTICAL POLE: 29mm diam., 400mm long
(1.14" diam. 15.74"L)

HORIZONTAL POLE: 25.4mm diam., 530mm long
(1.0" diam., 20.87"L)

MODEL S-4200: UNIVERSAL STAND WITH POST FOR F HOLDER

BASE: 265mm(W) x 265mm(D) x 25.4mm(H)
(10.43" x 10.43" x 1.0")

VERTICAL POLE: 29mm diam., 610mm long
(1.14" diam., 24.0"L)

HORIZONTAL POLE: 25.4mm diam., 530mm long
(1.0" diam., 20.86"L)

MODEL S-4300: UNIVERSAL STAND WITH ADAPTER FOR FS HOLDER

BASE: 265mm(W) x 265mm(D) x 25.4mm(H)
(10.43" x 10.43" x 1.0")

VERTICAL POLE: 29mm dia., 400mm long
(1.14" diam. 15.74"L)



HORIZONTAL POLE: 25.4mm diam., 535mm long
(1.0" diam., 21.06"L)

MODEL S-4400: UNIVERSAL STAND WITH ADAPTER FOR FS HOLDER

BASE: 265mm(W) x 265mm(D) x 25.4mm(H)
(10.43" x 10.43" x 1.0")

VERTICAL POLE: 29mm diam., 610mm long
(1.14" diam. 24.0"L)

HORIZONTAL POLE: 25.4mm diam., 535mm long
(1.0" diam., 21.06"L)



Model No.	Holder to Use	Vertical Pole	Horizontal Pole	Base
S-4100	F Holder	400mm (15.74")	530mm (20.86")	265mm x 265mm x 25.4mm (10.43" x 10.43" x 1.0")
S-4200	F Holder	610mm (24.0")	530mm (20.86")	
S-4300	FS Holder	400mm (15.74")	535mm (21.06")	
S-4400	FS Holder	610mm (24.0")	535mm (21.06")	

Material & Finish:

Base and Vertical/Horizontal Poles made of steel. Base finish in baked painting and Vertical / Horizontal poles rust-proofed in hard chrome. Slider block made of cast aluminum alloy in baked painting.



MEIJI TECHNO AMERICA

5895 Rue Ferrari, San Jose, CA 95138, USA

Tel : 408-226-3454 Fax : 408-226-0900

Toll Free : 1-800-832-0060

E-mail : info@meijitechno.com

http://www.meijitechno.com



MEIJI TECHNO CO., LTD.

322-1, Chikumazawa Miyoshi-machi, Iruma-gun
Saitama 354-0043, Japan

E-mail : meiji@meijitechno.co.jp

Phone : (0)49-259-0111 Fax : (0)49-259-0113

Cable : MEIJITECHNO OI