

# KBE-KBL Stand


Why settle for a small and inefficient work area under your microscope? MEIJI's new laminate base stand offer a clear, flat, 380mm x 330mm (15" x 13-1/4") work surface. The large work surface makes it easier to accommodate a multitude of objects. A variety of stereo microscope bodies and illuminators can be selected to fit your particular needs.



MODEL: EMZ-5/MA517/KBE/MA502

## KBE Stand

The focus block and holder are connected to a rigid extendable arm which is mounted to the flat laminate base. The picture shows the EMZ-5 Zoom stereo body with 0.5X supplementary lens mounted on the KBE Stand.



MODEL: EMT-2/MA264/KBL/MA502

## KBL Stand

The focus block and holder are mounted to a vertical 210mm (8") pillar which is attached to the flat laminate base. Optional 100mm (4") pillar extensions (MA551) can provide additional height and clearance.

The picture shows the EMT-2 Stereo body with the MA264 illuminator mounted on the KBL Stand.

Agent:

### MEIJI TECHNO CO.,LTD.

322-1,Chikumazawa Miyoshi-machi, Iruma-gun  
Saitama 354-0043, Japan  
E-mail : meiji@meijitechno.co.jp  
Phone : (0)49-259-0111 Fax : (0)49-259-0113



### MEIJI TECHNO AMERICA

5895 Rue Ferrari, San Jose, CA 95138, USA  
Tel : 408-226-3454 Fax : 408-226-0900  
Toll Free : 1-800-832-0060  
E-mail : info@meijitechno.com  
<http://www.meijitechno.com>

