


# New Product

# UL Stand




## UL STAND

Large universal stand featuring a pillar and crossarm mounted on a laminate board of 580mm x 600mm (22-1/3" x 19-1/4") with a rotatable stand for a light source. The pillar height is 260mm and the crossarm length is adjustable to 380mm (15").



(UL Stand does not include F Holder as shown.)



The above picture features MEIJI's SMD Inspection System composed of the following items:

- |                 |   |
|-----------------|---|
| <b>EMZ-5TR</b>  | Trinocular Zoom Stereo body   |
| <b>MA503</b>    | SWF15X Super wide field eyepieces (paired)  |
| <b>MA571</b>    | SMD Attachment (includes built-in 0.3X lens for oblique observation and 0.44X for vertical observation) |
| <b>F</b>        | Holder and focusing block   |
| <b>UL</b>       | Large universal stand with rotatable stand for light source   |
| <b>MA574</b>    | Sliding stage with 300mm x 200mm (12" x 8") X-Y travel  |
| <b>FL181</b>    | Dual arm fiber illuminator system   |
| <b>CK3600</b>   | Color TV camera, 1/2" CCD, 420 Horizontal TV lines, NTSC system   |
| <b>MA151/5N</b> | "C" mount adapter<br>14" Color TV Monitor   |