



MA507

1.5X Auxiliary Lens

Working Distance (49mm)



MA507 Auxiliary Lens 1.5X W.D. 49mm for selected EMZ bodies

KEY FEATURES

- Auxiliary Lens - 1.5X
- Working Distance: 49mm
- Fits: EMZ-1, EMZ-2, EMZ-2TR, EMZ-5, EMZ-5D, EMZ-5TR, EMZ-5TRD, EMZ-13, EMZ-13TR
- Warranty: LIMITED LIFETIME WARRANTY
- Manufacturer: MADE IN JAPAN



Fits: EMZ-5, EMZ-5D, EMZ-5TR, EMZ-5TRD, EMZ-13, EMZ-13TR bodies

SPECIFICATIONS

The MA507 is an auxiliary lens with 1.5X magnification and has a working distance of 49mm. Compatible with EMZ-1, EMZ-2, EMZ-2TR, EMZ-5, EMZ-5D, EMZ-5TR, EMZ-5TRD, EMZ-13, EMZ-13TR.

NOTE: Additional items are needed to create a full microscope system. The following components are normally ordered:

- Stereo Microscope Body (EMZ, EMT or EMF)
- Eyepieces (Requires 2 of identical magnification)
- Stand (pole stands, rigid arm stands, boom stands, heavy duty boom stands, wide surface stands or articulated arm stands)
- Illumination Source (Fiber optic halogen, fiber optic LED, fluorescent ring or LED illumination)
- Auxiliary Lens (Optional)

ADDITIONAL AUXILIARY LENS (OPTIONAL)

- **MA517:** 0.5X Auxiliary Lens W.D. 150mm (5.9")
- **MA526:** 0.75X Auxiliary Lens W.D. 97mm (3.8")

OPTIONAL EMZ EYEPIECES

- **MA501:** Super Widefield Eyepiece 5X, Field No. 26.0 (paired) (no reticle mount)
- **MA502:** Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- **MA718:** Super Widefield Eyepiece 7.5X, Field No. 26.0 (paired) (no reticle mount)
- **MA502:** Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- **MA520:** Super Widefield Eyepiece 12.5X, Field No. 20.0 (paired) (25mm reticle mount)
- **MA503:** Super Widefield Eyepiece 15X, Field No. 15.4 (paired) (25mm reticle mount)
- **MA535:** Super Widefield, High Eyepoint HSWF 15X Eyepiece (paired) (25mm reticle mount)
- **MA504:** Super Widefield Eyepiece 20X, Field No. 11.5 (paired) (25mm reticle mount)
- **MA521:** Super Widefield Eyepiece 30X, Field No. 7.7 (paired) (no reticle mount)
- **MA519:** Super Widefield Focusing Eyepiece 10X, Field No. 21 (each) (25mm reticle mount)

PACKAGING INFORMATION

- Ships in 8" x 8" x 5" box
- Weighs 0.5 lbs