



MA515

Adapter Used to Attach Auxiliary Lens for EMT Models



MA515 adapter used to attach MA518, MA527, MA513 and MA514 to selected EMT models

SPECIFICATIONS

The MA515 is an adapter used to attach MA518, MA527, MA513 and MA514 to EMT-2, EMTR-2, EMT-3, EMTR-3, EMT-4, and EMTR-4.

NOTE: Additional items are needed to create a full microscope system. The following components are normally ordered:

- Stereo Microscope Body (EMZ, EMT or EMF)
- Eyepieces (Requires 2 of identical magnification)
- Stand (pole stands, rigid arm stands, boom stands, heavy duty boom stands, wide surface stands or articulated arm stands)
- Illumination Source (Fiber optic halogen, fiber optic LED, fluorescent ring or LED illumination)
- Auxiliary Lens (Optional)

ADDITIONAL AUXILIARY LENS (OPTIONAL)

- **MA518:** Auxiliary Lens 0.5X, W.D. 108mm (4.3")
- **MA513:** Auxiliary Lens 1.5X, W.D. 35mm (1.4")
- **MA527:** Auxiliary Lens 0.75X, W.D. 72mm (2.8")
- **MA515:** Adapter used to attach MA518, MA527, MA513 to EMT-2

OPTIONAL EMZ/EMT EYEPIECES

- **MA501:** Super Widefield Eyepiece 5X, Field No. 26.0 (paired) (no reticle mount)
- **MA502:** Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- **MA718:** Super Widefield Eyepiece 7.5X, Field No. 26.0 (paired) (no reticle mount)
- **MA502:** Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- **MA520:** Super Widefield Eyepiece 12.5X, Field No. 20.0 (paired) (25mm reticle mount)
- **MA503:** Super Widefield Eyepiece 15X, Field No. 15.4 (paired) (25mm reticle mount)
- **MA535:** Super Widefield, High Eyepoint HSWF 15X Eyepiece (paired) (25mm reticle mount)
- **MA504:** Super Widefield Eyepiece 20X, Field No. 11.5 (paired) (25mm reticle mount)
- **MA521:** Super Widefield Eyepiece 30X, Field No. 7.7 (paired) (no reticle mount)
- **MA519:** Super Widefield Focusing Eyepiece 10X, Field No. 21 (each) (25mm reticle mount)

PACKAGING INFORMATION

- Ships in 8" x 8" x 5" box
- Weighs 0.3 lbs

KEY FEATURES

- Adapter for MA518, MA527, MA513 and MA514 to EMT-2, EMTR-2, EMT-3, EMTR-3, EMT-4, EMTR-4
- Warranty: LIMITED LIFETIME WARRANTY
- Manufacturer: MADE IN JAPAN



Fits: EMT-2, EMTR-2, EMT-3, EMTR-3, EMT-4 and EMTR-4