



# MA526

0.75X Auxiliary Lens

Working Distance (97mm)



MA526 Auxiliary Lens 0.75X W.D. 97mm for selected EMZ bodies

## KEY FEATURES

- Auxiliary Lens - 0.75X
- Working Distance: 97mm
- Fits: EMZ-1, EMZ-2, EMZ-2TR, EMZ-5, EMZ-5D, EMZ-5TR, EMZ-5TRD, EMZ-13, EMZ-13TR
- Warranty: LIMITED LIFETIME WARRANTY
- Manufacturer: MADE IN JAPAN



Fits: EMZ-5, EMZ-5D, EMZ-5TR, EMZ-5TRD, EMZ-13, EMZ-13TR bodies

## SPECIFICATIONS

The MA526 is an auxiliary lens with 0.75X magnification and has a working distance of 97mm. Compatible with EMZ-1, EMZ-2, EMZ-2TR, EMZ-5, EMZ-5D, EMZ-5TR, EMZ-5TRD, EMZ-13, EMZ-13TR.

**NOTE:** Additional items are needed to create a full microscope system. The following components are normally ordered:

- Stereo Microscope Body (EMZ, EMT or EMF)
- Eyepieces (Requires 2 of identical magnification)
- Stand (pole stands, rigid arm stands, boom stands, heavy duty boom stands, wide surface stands or articulated arm stands)
- Illumination Source (Fiber optic halogen, fiber optic LED, fluorescent ring or LED illumination)
- Auxiliary Lens (Optional)

## ADDITIONAL AUXILIARY LENS (OPTIONAL)

- MA507: 1.5X Auxiliary Lens W.D. 49mm (1.9")
- MA517: 0.5X Auxiliary Lens W.D. 150mm (5.9")

## OPTIONAL EMZ EYEPIECES

- MA501: Super Widefield Eyepiece 5X, Field No. 26.0 (paired) (no reticle mount)
- MA502: Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- MA718: Super Widefield Eyepiece 7.5X, Field No. 26.0 (paired) (no reticle mount)
- MA502: Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- MA520: Super Widefield Eyepiece 12.5X, Field No. 20.0 (paired) (25mm reticle mount)
- MA503: Super Widefield Eyepiece 15X, Field No. 15.4 (paired) (25mm reticle mount)
- MA535: Super Widefield, High Eyepoint HSWF 15X Eyepiece (paired) (25mm reticle mount)
- MA504: Super Widefield Eyepiece 20X, Field No. 11.5 (paired) (25mm reticle mount)
- MA521: Super Widefield Eyepiece 30X, Field No. 7.7 (paired) (no reticle mount)
- MA519: Super Widefield Focusing Eyepiece 10X, Field No. 21 (each) (25mm reticle mount)

## PACKAGING INFORMATION

- Ships in 8" x 8" x 5" box
- Weighs 0.5 lbs