



# MA527

0.75X Auxiliary Lens

Working Distance (72mm)



MA527 Auxiliary Lens 0.75X W.D. 72mm for EMT-2, EMTR-2, EMT-3, EMTR-3, EMT-4 and EMTR-4

## KEY FEATURES

- Auxiliary Lens - 0.75X
- Working Distance: 72mm
- Fits: EMT-2, EMTR-2, EMT-3, EMTR-3, EMT-4 and EMTR-4
- Warranty: LIMITED LIFETIME WARRANTY
- Manufacturer: MADE IN JAPAN



Fits: EMT-2, EMTR-2, EMT-3, EMTR-3, EMT-4 and EMTR-4

## SPECIFICATIONS

The MA527 is an auxiliary lens with 0.75X magnification and has a working distance of 72mm. Compatible with EMT-2, EMTR-2, EMT-3, EMTR-3, EMT-4, and EMTR-4.

**NOTE:** (MA515 adapter is needed to use MA527 on EMT-2)

**NOTE:** Additional items are needed to create a full microscope system. The following components are normally ordered:

- Stereo Microscope Body (EMZ, EMT or EMF)
- Eyepieces (Requires 2 of identical magnification)
- Stand (pole stands, rigid arm stands, boom stands, heavy duty boom stands, wide surface stands or articulated arm stands)
- Illumination Source (Fiber optic halogen, fiber optic LED, fluorescent ring or LED illumination)
- Auxiliary Lens (Optional)

## ADDITIONAL AUXILIARY LENS (OPTIONAL)

- **MA518:** Auxiliary Lens 0.5X, W.D. 108mm (4.3")
- **MA513:** Auxiliary Lens 1.5X, W.D. 35mm (1.4")
- **MA515:** Adapter used to attach MA518, MA527, MA513 to EMT-2

## OPTIONAL EMZ/EMT EYEPIECES

- **MA501:** Super Widefield Eyepiece 5X, Field No. 26.0 (paired) (no reticle mount)
- **MA502:** Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- **MA718:** Super Widefield Eyepiece 7.5X, Field No. 26.0 (paired) (no reticle mount)
- **MA502:** Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- **MA520:** Super Widefield Eyepiece 12.5X, Field No. 20.0 (paired) (25mm reticle mount)
- **MA503:** Super Widefield Eyepiece 15X, Field No. 15.4 (paired) (25mm reticle mount)
- **MA535:** Super Widefield, High Eyepoint HSWF 15X Eyepiece (paired) (25mm reticle mount)
- **MA504:** Super Widefield Eyepiece 20X, Field No. 11.5 (paired) (25mm reticle mount)
- **MA521:** Super Widefield Eyepiece 30X, Field No. 7.7 (paired) (no reticle mount)
- **MA519:** Super Widefield Focusing Eyepiece 10X, Field No. 21 (each) (25mm reticle mount)

## PACKAGING INFORMATION

- Ships in 8" x 8" x 5" box
- Weighs 0.5 lbs