



MA652

0.7X Auxiliary Lens

Working Distance (233mm)



MA652 Auxiliary Lens 0.7X W.D. 233mm for EMZ-12, EMZ-12D, EMZ-12TR, and EMZ-12TRD

KEY FEATURES

- Auxiliary Lens - 0.7X
- Working Distance: 233mm
- Fits: EMZ-12 Series
- EMZ-12 Threads: Thread M63.5, Pitch 0.75
- Warranty: LIMITED LIFETIME WARRANTY
- Manufacturer: MADE IN JAPAN



Fits: EMZ-12, EMZ-12D, EMZ-12TR, and EMZ-12TRD

SPECIFICATIONS

The MA652 is an auxiliary lens with 0.7X magnification and has a working distance of 233mm. Compatible with EMZ-12, EMZ-12D, EMZ-12TR and EMZ-12TRD.

NOTE: Additional items are needed to create a full microscope system. The following components are normally ordered:

- Stereo Microscope Body (EMZ, EMT or EMF)
- Eyepieces (Requires 2 of identical magnification)
- Stand (pole stands, rigid arm stands, boom stands, heavy duty boom stands, wide surface stands or articulated arm stands)
- Illumination Source (Fiber optic halogen, fiber optic LED, fluorescent ring or LED illumination)
- Auxiliary Lens (Optional)

PROTECTIVE GLASS (OPTIONAL)

- MA659: Non-shattering protective glass

OPTIONAL EMZ EYEPIECES

- MA501: Super Widefield Eyepiece 5X, Field No. 26.0 (paired) (no reticle mount)
- MA502: Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- MA718: Super Widefield Eyepiece 7.5X, Field No. 26.0 (paired) (no reticle mount)
- MA502: Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- MA520: Super Widefield Eyepiece 12.5X, Field No. 20.0 (paired) (25mm reticle mount)
- MA503: Super Widefield Eyepiece 15X, Field No. 15.4 (paired) (25mm reticle mount)
- MA535: Super Widefield, High Eyepoint HSWF 15X Eyepiece (paired) (25mm reticle mount)
- MA504: Super Widefield Eyepiece 20X, Field No. 11.5 (paired) (25mm reticle mount)
- MA521: Super Widefield Eyepiece 30X, Field No. 7.7 (paired) (no reticle mount)
- MA519: Super Widefield Focusing Eyepiece 10X, Field No. 21 (each) (25mm reticle mount)

PACKAGING INFORMATION

- Ships in 8" x 8" x 5" box
- Weighs 0.5 lbs