



# MA794

## 0.5X Auxiliary Lens


Working Distance (174mm)



MA794 Auxiliary Lens 0.5X W.D. 174mm for EMZ-8U, EMZ-8UD, EMZ-8UH, EMZ-8TR, EMZ-8TRD, EMZ-8TRH, and EMZ-8TRU

## KEY FEATURES

- Auxiliary Lens - 0.5X
- Working Distance: 174mm
- Fits: EMZ-8U and EMZ-8TR Series
- Warranty: LIMITED LIFETIME WARRANTY
- Manufacturer: MADE IN JAPAN



Fits: EMZ-8U and EMZ-8TR bodies

## SPECIFICATIONS

The MA794 is an auxiliary lens with 0.5X magnification and has a working distance of 174mm. Compatible with EMZ-8U, EMZ-8UD, EMZ-8UH, EMZ-8TR, EMZ-8TRD, EMZ-8TRH, and EMZ-8TRU.

**NOTE:** Additional items are needed to create a full microscope system. The following components are normally ordered:

- Stereo Microscope Body (EMZ, EMT or EMF)
- Eyepieces (Requires 2 of identical magnification)
- Stand (pole stands, rigid arm stands, boom stands, heavy duty boom stands, wide surface stands or articulated arm stands)
- Illumination Source (Fiber optic halogen, fiber optic LED, fluorescent ring or LED illumination)
- Auxiliary Lens (Optional)

## ADDITIONAL AUXILIARY LENS (OPTIONAL)

- MA792: 0.35X Auxiliary Lens W.D. 250mm (9.8")
- MA793: 0.44X Auxiliary Lens W.D. 184mm (7.2")
- MA795: 1.5X Auxiliary Lens W.D. 57mm (2.2")

## OPTIONAL EMZ EYEPIECES

- MA501: Super Widefield Eyepiece 5X, Field No. 26.0 (paired) (no reticle mount)
- MA502: Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- MA718: Super Widefield Eyepiece 7.5X, Field No. 26.0 (paired) (no reticle mount)
- MA502: Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- MA520: Super Widefield Eyepiece 12.5X, Field No. 20.0 (paired) (25mm reticle mount)
- MA503: Super Widefield Eyepiece 15X, Field No. 15.4 (paired) (25mm reticle mount)
- MA535: Super Widefield, High Eyepoint HSWF 15X Eyepiece (paired) (25mm reticle mount)
- MA504: Super Widefield Eyepiece 20X, Field No. 11.5 (paired) (25mm reticle mount)
- MA521: Super Widefield Eyepiece 30X, Field No. 7.7 (paired) (no reticle mount)
- MA519: Super Widefield Focusing Eyepiece 10X, Field No. 21 (each) (25mm reticle mount)

## PACKAGING INFORMATION

- Ships in 8" x 8" x 5" box
- Weighs 0.5 lbs