



# EMZ-10H

Binocular Greenough Design,  
Stereo Zoom Body (110mm W.D)  
High Eyepoint



EMZ-10H (0.7X–4.5X) Binocular Greenough Design, Stereo Zoom Body,  
Working Distance 4.3" (110mm)

## KEY FEATURES

- Zoom Ratio: 6.5:1
- Zoom Range: 0.7X–4.5X
- Magnification: 7X–45X (10X eyepieces)
- Field of View: 32mm–5.1mm
- Working distance: 110mm
- Extended Mag Range: 2.1X–270X (With auxiliary lenses and eyepieces)
- Mounting: Fits 84.5mm Diameter Focus Holders
- Warranty: LIMITED LIFETIME WARRANTY
- Manufacturer: MADE IN JAPAN

## SPECIFICATIONS

The EMZ-10H is a new and improved high eyepoint model with an eye relief of 20.9mm designed for microscopists and for people who wear glasses. It is specially designed to give the user more clearance between their eyes and the microscope. The high eyepoint models also uses a newly designed eyepiece with a wider designed lens to compliment the high eyepoint models (MA522 10X, MA818 15X, MA819 20X) Optional rubber eye guards (MA906) are also available.

The EMZ-10H body has a long 110mm working length, an increase focal depth and a wide field number which makes it perfect for all applications ranging from industry and biological fields to education and training. Despite these advanced capabilities, the microscope operation is simple. It features Greenough optics for enhance image flatness and contrast. As a result, views of three dimensional surfaces are sharp with excellent depth characteristic. The EMZ-10H offers high resolution observation in precision-oriented semi-conductor, popular for mounting onto bonders and probers, and offers superior performance at low cost. In addition to its compact design and versatility, a wide variety of auxiliary lenses and stands are available to allow the user to customize their desired working distance to adapt to almost any application. Suitable for eyeglass users.

**NOTE:** Additional items are needed to create a full microscope system. The following components are normally ordered:

- Stereo Microscope Body (Zoom, Turret or Fixed)
- Eyepieces (Requires 2 of identical magnification)
- Stand (pole stands, rigid arm stands, boom stands, heavy duty boom stands, wide surface stands or articulated arm stands)
- Illumination Source (Fiber optic halogen, fiber optic LED, fluorescent ring or LED illumination)
- Auxiliary Lens (Optional)

### Included with Microscope:

- MA906: Rubber eyeshields (paired)
- MA703: Vinyl dust cover for EM series

### Optional Eyepieces Available:

- MA522: Super Widefield Eyepiece High Eye Point 10X (paired)
- MA818/2: Super Widefield Eyepiece High Eye Point 15X (paired)
- MA819/2: Super Widefield Eyepiece High Eye Point 20X (paired)

### Optional Auxiliary Lens:

- MA558: Auxiliary Lens 0.4X, 247mm (9.7") W.D
- MA546: Auxiliary Lens 0.5X, 194mm (7.6") W.D
- MA547: Auxiliary Lens 0.75X, 127mm (5.0") W.D
- MA548: Auxiliary Lens 1.5X, 64mm (2.5") W.D

## PACKAGING INFORMATION

- Ships in 14" x 8" x 9" box
- Weighs 7 lbs