



F/BLACK Coarse

Focus Block/ Holder



F/BLACK coarse focus block has a 84.2mm diameter opening that accepts all EM Series stereo bodies. A donut adapter (MA1095) would be needed for stereo heads that are 76mm diameter.

KEY FEATURES

- Coarse Focus Mechanism
- Rack and Pinion Mount
- 55mm Travel
- 20mm Diameter
- 84.2mm diameter opening
- Color: Black
- Fits: All EM Series stereo bodies
- Warranty: LIMITED LIFETIME WARRANTY
- Manufacturer: MADE IN JAPAN

SPECIFICATIONS

F/BLACK Coarse Focus Block/Holder is a focusing holder with rack and pinion mount that allows smooth, precise focusing over 55mm of travel. Coarse focus mechanism incorporates a slip clutch and tension adjustment. Fits 20mm diameter, 71mm long mounting post. The 84.2mm diameter opening accepts all EM Series stereo bodies.

Focusing holders hold the microscope body and contain the mechanism that raises or lowers the body to focus on the specimen. Complete stands include their own focusing holders. But for other kinds of tasks, however, it may be necessary to mount the microscope body on a boom arm stand or directly onto a machine or some other surface.

NOTE: Additional items are needed to create a full microscope system. The following components are normally ordered:

- Stereo Microscope Body (EMZ, EMT or EMF)
- Eyepieces (Requires 2 of identical magnification)
- Stand (pole stands, rigid arm stands, boom stands, heavy duty boom stands, wide surface stands or articulated arm stands)
- Illumination Source (Fiber optic halogen, fiber optic LED, fluorescent ring or LED illumination)
- Auxiliary Lens (Optional)

Optional Focus Blocks Available

- **F:** Focusing slide block with coarse focus (Fits 20mm mounting post)
- **FC:** Focus block with coarse and fine focus (Fits 20mm mounting post, see FC11 below)
- **FC11:** Mounting pin extension (The FC11 is required when mounting the FC focus block on the 71mm long mounting pin of the S-4000 series stands) (FD and FCD Fits 20mm mounting post) (With 1/4"-20 mounting thread)
- **FD:** Focus block with coarse focus for mounting a CCD or digital camera (With 1/4"-20 mounting thread) (Fits 20mm mounting post)
- **FCD:** Focus block with coarse and fine focus for mounting a CCD or digital camera (With 1/4"-20 mounting thread) (Fits 20mm mounting post)

PACKAGING INFORMATION

- Ships in 12" x 9" x 7" box
- Weighs 4 lbs



F focus holder holding EMZ-5 and EMZ-13TRH bodies