



MA978

Rotatable 1/4 Wave Plate



Top




Bottom

SPECIFICATIONS

Rotatable 1/4 wave plate for used for EMZ body.

NOTE: Additional items are needed to create a full microscope system. The following components are normally ordered:

- Stereo Microscope Body (Zoom, Turret or Fixed)
- Eyepieces (Requires 2 of identical magnification)
- Stand (Pole stands, rigid arm stands, boom stands, heavy duty boom stands, wide surface stands or articulated arm stands)
- Illumination Source (Fiber optic halogen, fiber optic LED, fluorescent ring or LED illumination)
- Auxiliary Lens (Optional)



EMZ-13VX with MA978

PACKAGING INFORMATION

- Ships in 8" x 8" x 5" box
- Weighs 1 lb

KEY FEATURES

- Used with oblique incident-light illumination situations
- Size of Whole Wave Plate: 6cm in Diameter
- Warranty: LIMITED LIFETIME WARRANTY
- Manufacturer: MADE IN JAPAN