


Meiji Techno TC Series Inverted Microscopes

Mechanical Stage and Plate Inserts Dimensional Information



MA378 - 1 x 3 Glass Slide & Hemocytometer Holder
MA384 - 35mm Petri Dish Holder
MA387 - 65mm Petri Dish Holder
MA864 - Terasaki Holder

The condenser tower can be unscrewed and re-attached if room is needed to place a fixture on the stage.

All specimen holders use the same cutout thickness with the thinner dimension being 0.88mm to 0.95mm.

All specimen holder plates are 3mm thick and have the same outer dimensions.

