



DIGITAL CAMERA ADAPTERS




CANON
Powershot G3
G5 Bayonet type

CANON
Powershot G1,G2


NIKON
Coolpix 800, 900, 900S,
950, 990, 995, 4500

NIKON
Coolpix 5000


MA151/30/41




MA151/30/43




MA151/30/50




MA151/30/70




MA151/40/40




MA151/40/42




MA151/30/51




MA151/30/71




MA151/45/40




MA151/45/42




MA151/40/50




MA151/40/70




MA151/50/40




MA151/50/42



MA151/45/50



MA151/45/70




NIKON
 Coolpix 5400


NIKON
 Coolpix 5700

NIKON
 Coolpix 4300, 885


MA151/30/53




MA151/30/55




MA151/30/57




MA151/30/59




MA151/40/52




MA151/40/54




MA151/40/56




MA151/40/58




MA151/45/52




MA151/45/54




MA151/45/56




MA151/45/58




MA151/50/52




MA151/50/54



MA151/50/56



MA151/50/58



DIMENSIONAL ILLUSTRATIONS

media
0, C-3000
0, C-3100
0, C-5050


Olympus Camedia
C-700, C-720, C-740, C-750

FUJI FINEPIX
S602, 4900Z, 6900Z, S5000


SONY
S70, S75, S85
CD300, CD400, CD500

0/60


MA151/30/63



MA151/30/81




MA151/30/91




0/61


MA151/40/62



MA151/40/80




MA151/40/90




0/60


MA151/45/62



MA151/45/80




MA151/45/90




45/60


MA151/50/62



MA151/50/80



MA151/50/90



	NEW	NEW		
DIGITAL CAMERAS	CANON Powershot G3 G5 Bayonet type	CANON Powershot G1,G2	NIKON Coolpix 800, 900, 900S, 950 990 995 4500	NIKON Coolpix 5000
MICROSCOPES				
Camera thread in Diameter	49mm	46.5mm	28mm	46mm
Fits a 23.2mm I.D./25.2mm O.D. Eye tube or Photo tube of:	Adapter Cat.No.	Adapter Cat.No.	Adapter Cat.No.	Adapter Cat.No.
ML2000, ML5000, ML7000, ML8000, ML9000 Series EMTR-1, EMTR-2, EMTR-3, EMTR-4 EMZ-2TR, EMZ-5TR, EMZ-5TRD, EMZ-8TR. EMZ-8TRD, EMZ-12TR, EMZ-12TRD, GEMZ-8TR, GEMZ-5TR, GM-8TR, GM-5TR MA751 Photo attachment of RZ TM450, TM460, TM490 and other TM400 Series microscope with teaching head	MA151/30/41	MA151/30/43	MA151/30/50	MA151/30/70
Fits a 23.2mm I.D./27.2mm O.D. Eye tube of:				
TM201, TM202, TM203, TM204, TM205, TM206 TM207, TM208, TM290, TM210, TM211, TM212, TM213, TM214, TM215. TM420, TM440, TM420F, TM422F, TM440F, TM470PF, TM480P			MA151/30/51	MA151/30/71
Fits "C" mount thread, 25.4mm diameter				
MS-35, MS-35D, MS-35DR, MS-40, MS-40D, MS-40DR MS-45, MS-45D, MS-45DR, MS-50, MS-50D, MS-50DR MS-15 VM-1V, ML-1VL, VM-1D, VM-2V, VM-2VL	MA151/40/40	MA151/40/42	MA151/40/50	MA151/40/70
Fits a 30.5mm I.D./34.0mm O.D. Eye tube of:				
EMF-1, EMF-2, EMT-1, EMT-2, EMT-3, EMT-4, EMZ-1, EMZ-2, EMZ-5, EMZ-6, EMZ-9, EMZ-10, EMZ-12, Z-7000, GEMZ-5, GEMT-2, GEMT-4, GM-2T, GM-3T, GM-4T, SK and SKT Series	MA151/45/40	MA151/45/42	MA151/45/50	MA151/45/70
Fits a 30.0mm I.D., 34.0mm O.D.				
Any microscope with an image plane at 10mm from top end of eye tube.	MA151/50/40	MA151/50/42		

COUPLERS

MA151/37

Coupler to convert the M28 x 0.75 thread of MA151/40/50 to a M37 x 0.75 internal thread
This is used in conjunction with the MA151/40/50.

MA151/43

Coupler to convert the M28 x 0.75 thread of MA151/40/50 to a M43 x 0.75 internal thread
This is used in conjunction with the MA151/40/50.

How to use

For photo tube of trinocular head, slip the adapter over the top of the tube and secure it with locking screw.

For eye tube of binocular head, insert the adapter into the eye tube until fully seated and secure the adapter with lock down screw.



DIGITAL CAMERA ADAPTERS

NEW	NEW	NEW		NEW		NEW
NIKON Coolpix 5400	NIKON Coolpix 5700	NIKON Coolpix 4300, 885	Olympus Camedia C-2000,C-2020,C-3000 C-3030,C-3040,C-3100 C-4040,C-4100,C-5050	Olympus Camedia C-700, C-720,C-740 C-750	FUJI FINEPIX S602, 4900Z 6900Z, S5000	SONY S70, S75, S85 CD300, CD400 CD500
45mm	53.4mm	36.5mm	41mm	45.6mm	55mm	45mm
Adapter Cat.No.	Adapter Cat.No.	Adapter Cat.No.	Adapter Cat.No.	Adapter Cat.No.	Adapter Cat.No.	Adapter Cat.No.
MA151/30/53	MA151/30/55	MA151/30/57	MA151/30/61	MA151/30/63	MA151/30/81	MA151/30/91
MA151/40/52	MA151/40/54	MA151/40/56	MA151/40/60	MA151/40/62	MA151/40/80	MA151/40/90
MA151/45/52	MA151/45/54	MA151/45/56	MA151/45/60	MA151/45/62	MA151/45/80	MA151/45/90
MA151/50/52	MA151/50/54	MA151/50/56		MA151/50/62	MA151/50/80	MA151/50/90



Agent :

MEIJI TECHNO AMERICA
 2186 Bering Drive, San Jose, CA., 95131, USA
 Phone : 408-428-9654 Fax : 408-428-0472
 Toll free : 800-832-0060
 E-mail : sales@meijitechno.com



MEIJI TECHNO CO.,LTD.
 322-1,Chikumazawa Miyoshi-machi
 Iruma-gun Saitama 354-0043, Japan
 Phone : (0)49-259-0111 Fax : (0)49-259-0113
 E-mail : meiji@meijitechno.co.jp

<http://www.meijitechno.com>
<http://www.meijitechno.co.jp>