


Meiji Techno Boom Mounted Optics Part Numbers & Specifications



System# vmboom

This compound optics boom setup is the most flexible system for high magnification, high-resolution imaging. The stand can accommodate any size specimen while the 150 watt fiber optic epi-illumination provides bright, through-the-objective lighting to any surface or specimen.

Part Numbers & Specifications

Part Numbers	Description
A	S-4100 Universal Boom Stand with mounting post, requires focus block, vertical pole: 400mm, horizontal: 530mm. Base: 265mm x 265mm x 25mm.
B	No Eyepieces No Eyepieces for this configuration
C	VM-035 Video Microscope, 0.35X factor, with vertical Cmount & iris
	FL150/95 Single light guide, flexible, 36" length, 7.9mm end + 15mm end
D	VM-PC/SP1 VM mount with coarse and fine focus block minus "P" Stand
E	FT190/115 150 watt fiber optic light source and power supply
F	CC2100 USB2 Digital Video Camera, 1.4MP, complete with software



**MEIJI
TECHNO**

Meiji Techno America 3010 Olcott Street Santa Clara, CA 95054
Toll Free USA: 1(800) 832-0060 (408) 970-4799 FAX: (408) 970-5054
Visit us online at www.meijitechno.com