



MA504

Super Widefield Eyepiece 20X (Paired)



MA504 Super Widefield 20X Eyepiece designed for EMZ, EMT, EMF models

KEY FEATURES

- 20X Eyepieces
- Field No. 11.5
- 30.5mm Outer Diameter
- Accepts 25mm Reticle
- Made in Japan
- Limited Lifetime Warranty



A pair of optional rubber eyeshields are available for EMZ, EMT and EMF eyepieces (MA600).

SPECIFICATIONS

MA504 (Paired)

The MA504 is a Super Widefield High Eyepoint Eyepiece 20X (pair) Field No. 11.5, with 30.5mm O.D. (Outer Diameter), with 25mm reticle mount. This 20X Super Widefield Eyepiece is designed to fit our Meiji Techno models EMZ (Zoom Series), EMT (Turret Series), EMF (Fixed Series), but not compatible with High Eyepoint Models.

NOTE: Additional items are needed to create a full microscope system. The following components are normally ordered:

- Stereo Microscope Body (EMZ, EMT or EMF)
- Eyepieces (Requires 2 of identical magnification)
- Stand (pole stands, rigid arm stands, boom stands, heavy duty boom stands, wide surface stands or articulated arm stands)
- Illumination Source (Fiber optic halogen, fiber optic LED, fluorescent ring or LED illumination)
- Auxiliary Lens (Optional)

Optional Eyepieces Available

- **MA501:** Super Widefield Eyepiece 5X, Field No. 26.0 (paired) (no reticle mount)
- **MA718:** Super Widefield Eyepiece 7.5X, Field No. 26.0 (paired) (no reticle mount)
- **MA502:** Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- **MA520:** Super Widefield Eyepiece 12.5X, Field No. 20.0 (paired) (25mm reticle mount)
- **MA503:** Super Widefield Eyepiece 15X, Field No. 15.4 (paired) (25mm reticle mount)
- **MA535:** Super Widefield, High Eyepoint HSWF 15X Eyepiece (paired) (25mm reticle mount)
- **MA521:** Super Widefield Eyepiece 30X, Field No. 7.7 (paired) (no reticle mount)
- **MA519:** Super Widefield Focusing Eyepiece 10X, Field No. 21 (each) (25mm reticle mount)
- **MA600:** Rubber Eyeshields for SWF eyepieces (paired) (Fits 30.5mm O.D eyepieces) (fits EMZ, EMT, and EMF models)

PACKAGING INFORMATION

Ships in 8" x 8" x 5" package
Weight: 0.3 lbs