



MA906

Eyeshields for Super Widefield High Eyepoint Eyepieces (Paired)



MA906 Eyeshield for Super Widefield High Eyepoint Eyepieces

KEY FEATURES

- Rubber Eyeguards
- For EMZ High Eyepoint Models
- Fits over 30.0 mm O.D eyepieces
- Designed for MA522, MA818, and MA819 eyepieces
- Color: Matte Black
- Warranty: LIMITED LIFETIME WARRANTY
- Manufacturer: MADE IN JAPAN



MA906 on EMZ-5TRH and EMZ-13H High Eyepoint Models

SPECIFICATIONS

MA906 are rubber eyeshields (paired) for SWF eyepieces (each) (Fits MA522, MA818, and MA819).

For High Eyepoint models only (EMZ-5H, EMZ-5TRH, EMZ-10H, EMZ-8UH, EMZ-8TRH, EMZ-13H, EMZ-13TRH)

NOTE: Additional items are needed to create a full microscope system. The following components are normally ordered:

- Stereo Microscope Body (EMZ, EMT or EMF)
- Eyepieces (Requires 2 of identical magnification)
- Stand (pole stands, rigid arm stands, boom stands, heavy duty boom stands, wide surface stands or articulated arm stands)
- Illumination Source (Fiber optic halogen, fiber optic LED, fluorescent ring or LED illumination)
- Auxiliary Lens (Optional)

EMZ High Eyepoint Eyepieces

- MA522: Super Widefield Eyepiece High Eye Point 10X (paired)
- MA818/2: Super Widefield Eyepiece High Eye Point 15X (paired)
- MA819/2: Super Widefield Eyepiece High Eye Point 20X (paired)

EMZ Eyepieces (Not High Eyepoint)

- MA501: Super Widefield Eyepiece 5X, Field No. 26.0 (paired) (no reticle mount)
- MA718: Super Widefield Eyepiece 7.5X, Field No. 26.0 (paired) (no reticle mount)
- MA502: Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- MA520: Super Widefield Eyepiece 12.5X, Field No. 20.0 (paired) (25mm reticle mount)
- MA503: Super Widefield Eyepiece 15X, Field No. 15.4 (paired) (25mm reticle mount)
- MA535: Super Widefield, High Eyepoint HSWF 15X Eyepiece (paired) (25mm reticle mount)
- MA504: Super Widefield Eyepiece 20X, Field No. 11.5 (paired) (25mm reticle mount)
- MA521: Super Widefield Eyepiece 30X, Field No. 7.7 (paired) (no reticle mount)
- MA519: Super Widefield Focusing Eyepiece 10X, Field No. 21 (each) (25mm reticle mount)
- MA600: Rubber Eyeshields for SWF eyepieces (paired) (Fits 30.5mm O.D eyepieces) (fits EMZ, EMT, and EMF models)

PACKAGING INFORMATION

- Ships in 8" x 8" x 5" package
- Weight: 0.3 lbs