



EMT-2

Turret Stereo Body (1X and 3X)

Working Distance 81mm



EMT-2 (1.0X and 3.0X) Turret Stereo Body Working Distance 81mm

KEY FEATURES

- Magnification: 10X, 30X (10X eyepieces)
- Field of View: 23mm, 7.6mm
- Working distance: 81mm
- Extended Mag Range: 5X – 180X
- Auxiliary Lenses: 0.5X (108mm W.D.), 0.75X (72mm W.D.), 1.5X (35mm W.D.)
- Mounting: Fits 84.5mm Diameter Focus Holders
- Warranty: LIMITED LIFETIME WARRANTY
- Manufacturer: MADE IN JAPAN



A pair of optional rubber eyeshields are available for EMZ, EMT and EMF eyepieces (MA600).

SPECIFICATIONS

The EMT-2 turret stereo microscope body permits magnification to varied by means of a objective turret changer. The EMT-2 microscope body enable one-touch magnification power switch over between 1X and 3X. Because the operator's hand is kept at a distance, it decreases the possibilities the contamination due to dust and particulates is minimize. Its long 81mm working distance and a wide field of view makes the EMT-2 a perfect turret microscope head that provides excellent maneuverability proving optimum for special applications such as gemology, horology, hair transplants, and dental labs.

NOTE: Additional items are needed to create a full microscope system. The following components are normally ordered:

- Stereo Microscope Body (Zoom, Turret or Fixed)
- Eyepieces (Requires 2 of identical magnification)
- Stand (Pole stands, rigid arm stands, boom stands, heavy duty boom stands, wide surface stands or articulated arm stands)
- Illumination Source (Fiber optic halogen, fiber optic LED, fluorescent ring or LED illumination)
- Auxiliary Lens (Optional)

Included with Microscope:

- MA600: Rubber Eyeshields for SWF eyepieces (paired) (Fits 30.5mm O.D eyepieces)
- MA703: Vinyl dust cover for EM series

Optional Eyepieces Available:

- MA501: Super Widefield Eyepiece 5X, Field No. 26.0 (paired) (no reticle mount)
- MA502: Super Widefield Eyepiece 10X, Field No. 23.0 (paired) (25mm reticle mount)
- MA520: Super Widefield Eyepiece 12.5X, Field No. 20.0 (paired) (25mm reticle mount)
- MA503: Super Widefield Eyepiece 15X, Field No. 15.4 (paired) (25mm reticle mount)
- MA504: Super Widefield Eyepiece 20X, Field No. 11.5 (paired) (25mm reticle mount)

Optional Auxiliary Lens:

- MA518: Auxiliary Lens 0.5X, W.D. 108mm (4.3")
- MA527: Auxiliary Lens 0.75X, W.D. 72mm (2.8")
- MA513: Auxiliary Lens 1.5X, W.D. 35mm (1.4")
- MA515: Adapter used to attach MA518, MA527, MA513 to EMT-2

PACKAGING INFORMATION

- Ships in 14" x 8" x 9" box
- Weighs 5 lbs